cctbx PDB handling tools

Ralf W. Grosse-Kunstleve^a and Paul D. Adams^{a,b}

^aLawrence Berkeley National Laboratory, Berkeley, CA 94720

^bDepartment of Bioengineering, University of California at Berkeley, Berkeley, CA 94720

Correspondence email: RWGrosse-Kunstleve@LBL.Gov

Introduction

The PDB format is the predominant working format for atomic parameters (coordinates, occupancies, displacement parameters, etc.) in macromolecular crystallography. Many small-molecule programs also support this format. The PDB format specifications are available at http://www.pdb.org/. Technically, the format is very simple, therefore a vast number of parsers exist in scientific packages. This article is about the PDB handling tools included in the Computational Crystallography Toolbox (cctbx, http://cctbx.sourceforge.net/), the open-source component of the PHENIX project (http://www.phenix-online.org/).

The evolution of the cctbx PDB handling tools has gone through three main stages spread out over several years. A simple parser implemented in Python has been available for a long time. In many cases Python's runtime performance is sufficient for interactive processing of PDB files, but can be limiting for large files, or for repeatedly traversing the entire PDB database. This has prompted us to implement a fast C++ parser that is described in Grosse-Kunstleve et al. (2006). However, initially the fast cctbx PDB handling tools only supported "read-only" access. Writing of PDB files was supported only at a very basic level. This shortcoming has been removed and the current cctbx version provides comprehensive tools for reading, manipulating, and writing PDB files. These tools are available from both Python and C++, under the iotbx.pdb module.

This article presents an overview of the main types in the iotbx.pdb module, considerations that lead to the design, and related important nomenclature. It is not a tutorial for using the iotbx.pdb facilities. For this, refer to http://cctbx.sourceforge.net/sbgrid2008/tutorial.html. See also http://cci.lbl.gov/hybrid36/ which describes iotbx.pdb facilities for handling very large models.

Real-world PDB files

The simplicity of the PDB format is only superficial and, in the general case, stops after the initial parsing level. The structure of the PDB file implies a *hierarchy* of objects. A first approximation is this hierarchical view is:

```
model
chain
residue
atom
```

This is only an approximation because of a feature that is easily overlooked at first: the "altloc" (official PDB nomenclature) column 17 of PDB ATOM records, specifying "alternative location" identifiers for atoms in alternative conformations. As it turned out, about 90% of the development time invested into iotbx.pdb was in some form related to alternative conformations. Our goal was to provide robust tools that work even for the most unusual (but valid) cases, since this is a vital characteristic of any automated system. The main difficulties encountered while pursuing this goal were:

- Chains with conformers that have different sequences
- Chains with duplicate resseq+icode (residue sequence number + insertion code)
- Conformers interleaved or separated

These difficulties are best illustrated with examples. An old version of PDB entry 2IZQ (from Aug 2008) includes a chain with conformers that have different sequences, the residue with resseq 11 in chain A, atoms 220 through 283:

HEADER	AN	TIB	IOTIC					26-JU	JL-06	2IZQ	
ATOM	220	N	ATRP	Α	11	20.498	12.832	34.558	0.50	6.03	N
			.ATRP	Α	11						
ATOM	243	НН	2ATRP	Α	11	15.522	9.077	38.323	0.50	10.40	Н
ATOM	244	N	CPHE	Α	11	20.226	13.044	34.556	0.15	6.35	N
			.CPHE	Α	11						
ATOM	254	CZ	CPHE	Α	11	16.789	9.396	34.594	0.15	10.98	С
ATOM	255	N	BTYR	Α	11	20.553	12.751	34.549	0.35	5.21	N
			.BTYR	Α	11						
ATOM	261	CD	1BTYR	Α	11	18.548	10.134		0.35	9.45	С
ATOM	262	НВ	2CPHE	Α	11	21.221	10.536	34.146	0.15	7.21	Н
ATOM	263	CD	2BTYR	Α	11	18.463	10.012	36.681	0.35	9.08	С
ATOM	264	HB	3СРНЕ	Α	11	21.198	10.093	35.647	0.15	7.21	Н
ATOM	265	CE	1BTYR	Α	11	17.195	9.960	34.223	0.35	10.76	С
ATOM	266	HD	1CPHE	Α	11	19.394	9.937	32.837	0.15	10.53	Н
ATOM	267	CE	2BTYR	Α	11	17.100	9.826	36.693	0.35	11.29	С
ATOM	268	HD	2CPHE	Α	11	18.873	10.410	36.828	0.15	9.24	Н
ATOM	269	CZ	BTYR	Α	11	16.546	9.812	35.432	0.35	11.90	C
ATOM	270	HE	1CPHE	Α	11	17.206	9.172	32.650	0.15	12.52	Н
ATOM	271	ОН	BTYR	Α	11	15.178	9.650	35.313	0.35	19.29	0
ATOM	272	HE	2CPHE	Α	11	16.661	9.708	36.588	0.15	11.13	H
ATOM	273	HZ	CPHE	Α	11	15.908	9.110	34.509	0.15	13.18	Н
ATOM	274	H	BTYR	Α	11	20.634	12.539	33.720	0.35	6.25	Н
			.BTYR	Α	11						
ATOM	282	НН	BTYR	Α	11	14.978	9.587	34.520	0.35	28.94	Н

The original atom numbering does not have gaps. Here we have omitted blocks of atoms with constant resname and resseq+icode to save space.

As of Jun 8 2010, there are 74 files with mixed residue names in the PDB, i.e. only about 0.1% of the files. However, these files are perfectly valid and a PDB processing library is suitable as a component of an automated system only if it handles them sensibly.

An old version of PDB entry 1ZEH (Aug 2008) includes a chain with consecutive duplicate resseq+icode, atoms 878 through 894:

HEADER	HO	RMO	ΝE		01-MAY-98 1ZEH						
HETATM	878	C1	ACRS	5	12.880	14.021	1.197	0.50	33.23	С	
HETATM	879	C1	BCRS	5	12.880	14.007	1.210	0.50	34.27	С	
			.ACRS	5							
HETATM	892	01	ACRS	5	11.973	14.116	2.233	0.50	34.24	0	
HETATM	893	01	BCRS	5	11.973	14.107	2.248	0.50	35.28	0	
HETATM	894	0	НОН	5	-0.924	19.122	-8.629	1.00	11.73	0	
HETATM	895	0	НОН	6	-19.752	11.918	3.524	1.00	13.44	0	
HETATM	896	0	HOH	7	-1.169	17.936	-6.103	1.00	12.89	0	

To a human inspecting the old 1ZEH entry, it is of course immediately obvious that the water with resseq 5 is not related to the previous residue with resseq 5. However, arriving at this conclusion with an automatic procedure is not entirely straightforward. The human brings in the knowledge that water atoms without hydrogen are not covalently connected to other atoms. This is very detailed, specialized knowledge. Introducing such heuristics into an automatic procedure is likely to lead to surprises in some situations and is best avoided, if possible.

In the PDB archive, alternative conformers of a residue always appear consecutively. However, as mentioned in the introduction, the PDB format is also the predominant working format. Some programs produce files with conformers separated in this way (this file was provided to us by a user):

```
1716
              N
                  ALEU
                          190
                                    28.628
                                               4.549
                                                       20.230
                                                                                        N
ATOM
        1717
               CA ALEU
                          190
                                    27.606
                                               5.007
                                                       19.274
                                                                0.70
                                                                       3.71
                                                                                        С
ATOM
        1718
                          190
                                                       18.800
                                                                0.70
                                                                                        С
               CB ALEU
                                    26.715
                                               3.852
                                                                       4.15
ATOM
        1719
               CG ALEU
                          190
                                    25.758
                                               4.277
                                                       17.672
                                                                0.70
                                                                       4.34
                                                                                        С
        1829
                          190
                                                       20.343
                                                                0.30
                                                                                        N
ATOM
              N
                  BLEU
                                    28.428
                                               4.746
                                                                      5.13
                                                                                        С
ATOM
        1830
               CA BLEU
                          190
                                    27.378
                                               5.229
                                                       19.418
                                                                0.30
                                                                       4.89
                                                                                        С
ATOM
        1831
               CB
                  BLEU
                          190
                                    26.539
                                               4.062
                                                       18.892
                                                                0.30
                                                                       4.88
ATOM
        1832
               CG
                  BLEU
                          190
                                    25.427
                                               4.359
                                                       17.878
                                                                0.30
                                                                      5.95
                                                                                        С
                                                                                        N
ATOM
        1724
              N
                  ATHR
                          191
                                    27.350
                                               7.274
                                                       20.124
                                                                0.70
                                                                       3.35
                                                                                        С
ATOM
        1725
               CA ATHR
                          191
                                    26.814
                                               8.243
                                                       21.048
                                                                0.70
                                                                       3.27
                                                                                        С
ATOM
        1726
               CB ATHR
                          191
                                    27.925
                                               9.229
                                                       21.468
                                                                0.70
                                                                       3.73
                                                                                        0
ATOM
        1727
               OG1ATHR
                          191
                                    28.519
                                               9.718
                                                       20.259
                                                                0.70
                                                                      5.22
ATOM
        1728
               CG2ATHR
                          191
                                    28.924
                                               8.567
                                                       22.345
                                                                0.70
                                                                                        С
                                                                       4.21
ATOM
        1729
              C
                  ATHR
                          191
                                    25.587
                                               8.983
                                                       20.559
                                                                0.70
                                                                       3.53
                                                                                        С
ATOM
        1730
              0
                  ATHR
                          191
                                    24.872
                                               9.566
                                                       21.383
                                                                0.70
                                                                       3.93
                                                                                        0
                          ... residues 191 through 203 not shown
ATOM
        1828
              0
                  AGLY
                          203
                                             14.861
                                                       23.401
                                                                0.70
                                                                      5.84
                                                                                        0
                                     8.948
ATOM
        1833
              CD1BLEU
                          190
                                    26.014
                                               4.711
                                                       16.521
                                                                0.30
                                                                       6.21
                                                                                        C
ATOM
        1835
              C
                  BLEU
                          190
                                    26.506
                                               6.219
                                                       20.135
                                                                0.30
                                                                       4.99
                                                                                        C
ATOM
        1836
              0
                  BLEU
                          190
                                    25.418
                                               5.939
                                                       20.669
                                                                0.30
                                                                      5.91
                                                                                        0
        1721
                          190
                                    24.674
                                                                0.70
                                                                      5.31
                                                                                        C
ATOM
              CD2ALEU
                                               3.225
                                                       17.536
        1722
                          190
                                               6.055
                                                       20.023
                                                                0.70
                                                                                        С
ATOM
              C
                  ALEU
                                    26.781
                                                                       3.36
                                                       20.563
                                                                0.70
        1723
              0
                          190
                                    25.693
                                               5.796
                                                                                        0
ATOM
                  ALEU
                                                                      3.68
              С
                          190
                                                       20.023
                                                                                        C
ATOM
        8722
                  DLEU
                                    26.781
                                               6.055
                                                                0.70
                                                                       3.36
ATOM
        8723
              0
                  DLEU
                          190
                                    25.693
                                               5.796
                                                       20.563
                                                                0.70
                                                                       3.68
                                                                                        0
                                                       20.023
ATOM
        9722
              C
                          190
                                                                0.70
                                                                                        C
                  CLEU
                                    26.781
                                               6.055
                                                                       3.36
ATOM
        9723
              0
                  CLEU
                          190
                                    25.693
                                               5.796
                                                       20.563
                                                                0.70
                                                                       3.68
```

In this file, conformers A and B of residue 190 appear consecutively, but conformers C and D appear only after conformers A and B of all residues 191 through 203. While this is not the most intuitive way of ordering the residues in a file, it is still considered valid because the original intention is clear. Since it was our goal to develop a comprehensive library suitable for automatically processing files produced by any popular program, we found it important to correctly handle non-consecutive conformers.

iotbx.pdb.hierarchy

When developing the procedure capable of handling the variety of real-world situations shown above, we strived to keep the underlying rule-set as simple as possible and to avoid highly specific heuristics (e.g. "water is never covalently bound"). Complex rules imply complex implementations, are difficult to explain and understand, tend to lead to surprises, and are therefore likely to be rejected by the community. With this and the real-world situations in mind, we arrived at the following *primary* organization of the PDB hierarchy:

```
Primary PDB hierarchy:
    model(s)
    id
    chain(s)
    id
    residue_group(s)
    resseq
    icode
    atom_group(s)
    resname
    altloc
    atom(s)
```

In this presentation the "(s)" indicates a list of objects of the given type. i.e. a hierarchy contains a list of models, each model has an "id" (a simple string) and holds a list of chains, etc.

Comparing with the "first approximation hierarchy" above, the residue type is replaced with two new types: residue_group and atom_group. These types had to be introduced to cover all the real-world cases shown above. The residue_group and atom_group types are new and unusual. Before we go into the details of these types, it will be helpful to consider the bigger picture by introducing the alternative *secondary* view of the PDB hierarchy:

```
Secondary view of PDB hierarchy:
  model(s)
   id
   chain(s)
    id
   conformer(s)
     altloc
   residue(s)
     resname
   resseq
   icode
   atom(s)
```

This organization is probably more intuitive at first. The first two levels (model, chain) are exactly the same as in the primary hierarchy. Each chain holds a list of conformer objects, which are characterized by the altloc character from column 17 in the PDB ATOM records. A conformer is understood to be a *complete copy* of a chain, but usually two conformers *share* some or even most atoms. A residue is characterized by a unique resname and the resseq+icode.

The secondary view of the hierarchy evolved in the context of generating geometry restraints for refinement (e.g. bond, angle, and dihedral restraints), where this organization is most useful. It is also the organization introduced in Grosse-Kunstleve et al. (2006), where it was actually the primary organization. While working with the conformer-residue organization, we found that it is difficult to manipulate a hierarchy (e.g. add or delete atoms) in obvious ways. Finally, while developing the automatic generation of constrained occupancy groups for alternative conformations, the conformer-residue organization proved to be unworkable. The main difficulty is that the relative order of residues with alternative conformations is lost in the conformer-residue organization; it is only given indirectly by interleaved residues with shared atoms -- if they exist. As convenient as the conformer-residue organization is for the generation of restraints, it is a hindrance for other purposes.

The names for the new types in the primary hierarchy were chosen not to collide with the secondary view. We could have used "conformer" and "residue" again, just reversed, but there would be the big surprise that one residue has different resnames. To send the signal "this is not what you usually think of as a residue", we decided to use residue_group as the type name. A residue group holds a list of atom_group objects. All atoms in an atom group have the same resname and altloc. Therefore "resname_altloc_group" would have been another plausible name, but we favored atom group since it is more concise and better conveys what is the main content.

Detection of residue groups and atom groups

When constructing the primary hierarchy given a PDB file, the processing algorithm has to detect models, chains, residue groups, atom groups and atoms. Most steps are fairly straightforward, but none of the steps is actually completely trivial. For example, what to do if TER or ENDMDL cards are missing? What if residue sequence numbers are not consecutive? The ribosome community widely uses segid instead of chain id (even though the segid column is officially deprecated by the PDB). What if a file

contains both chain id and segid? Documenting all our answers in full detail is beyond the scope of this article (and would be more distracting than helpful anyway because the source code is openly available). Therefore we concentrate on the most important rules for the detection of residue groups and atom groups.

Then central conflict we had to resolve was:

- In order to handle non-consecutive conformers, we have to use the resseq+icode to find and group related residues.
- However, we cannot use the resseq+icode alone as a guide, because we also want to handle chains with duplicate resseq+icode (consecutive or non-consecutive).

This lead to a two stage procedure. In the first stage:

• A residue group is given by a block of consecutive ATOM or HETATM records with identical resseq+icode columns (SIGATM, ANISOU, and SIGUIJ records may be interleaved), e.g.

ATOM	234	H	ATRP	Α	11	20.540	12.567	33.741	0.50	7.24	H
MOTA	235	HA	ATRP	Α	11	20.771	12.306	36.485	0.50	6.28	H
ATOM	244	N	CPHE	Α	11	20.226	13.044	34.556	0.15	6.35	N
ATOM	245	CA	CPHE	Α	11	20.950	12.135	35.430	0.15	5.92	C

However, there is an important pre-condition:

• Unless a sub-block of ATOM records with identical resname+resseq+icode columns contains a main-conformer atom (almost "blank altloc"), e.g.

HEADER	Al	NTIBI	OTIC	RESIST	TANCE	07-MAY-97 1AJQ						
CRYST1	52	.120	65.	080	76.300 100.2	0 111.44	105.81	P 1		1		
HETATM	6097	CA	CA	1	5.676	34.115	52.446	1.00	18.50		CA	
HETATM	6100	C2	SPA	1	11.860	36.159	33.853	1.00	14.30		С	
HETATM	6107	C6	SPA	1	13.085	36.522	34.644	1.00	17.34		С	

In this case the sub-block is assigned to a separate residue group.

Within each residue group, all atoms are grouped by altloc+resname and assigned to atom groups. The order of the atoms in a residue group does not affect the assignment to atom groups.

After all atom records are assigned to residue groups and atom groups,

- residue groups with identical resseq+icode
- that do not contain main-conformer atoms

are merged in the second processing stage. While two residue groups are merged, atom groups with identical altloc+resname from the two sources are also merged.

Note that the grouping steps in this process may change the order of the atoms. However, our implementation preserves the relative order of the atoms as much as possible. I.e. in the hierarchy, resseq+icode appear in the original "first seen" order, and similarly for altloc+resname within a residue group.

Construction of secondary view of hierarchy

The secondary view of the hierarchy is constructed trivially from the primary hierarchy objects. For each chain, the complete list of altloc characters is determined in a first pass. In a second pass, a conformer object is created for each altloc, and a second loop over the chain assigns the atoms to each conformer. Main-conformer atoms are assigned to all conformers, atoms in alternative

conformations only to the corresponding conformer. Because of the difficulties alluded to earlier, the conformer and residue objects in the secondary view are "read-only". I.e. all manipulations such as addition or removal of residues, have to be performed on the primary hierarchy. Re-constructing the secondary view after the primary hierarchy has been changed is very fast (fractions of seconds even for the largest files, e.g. 0.22 s for PDB entry 1HTQ with almost one million atoms).

The iotbx/examples/pdb_hierarchy.py script shows how to construct the primary hierarchy from a PDB file, and how to obtain the secondary view.

Nomenclature related to alternative conformations

The iotbx.pdb module uses the following nomenclature to describe various aspects of alternative conformations:

- *Main conf. atom* : an atom with
 - o a blank altloc character
 - o and no other atom with the same name+resname (but different altloc) in the residue group. This second condition is needed for cases like this:

HEADER	HY	DROLASE					12-MA	AR-04	1S07	
ATOM	2460	CG1 VAL	Α	325	-23.284	97.713	15.815	0.66	21.74	С
ATOM	2461	CG1AVAL	Α	325	-23.010	97.616	18.295	0.66	22.88	C
ATOM	2462	CG2BVAL	Α	325	-24.819	96.373	17.146	0.66	22.57	C

- Alt. conf. atom: not main conf.
- Conformer: a complete chain with main conf. atoms and alt. conf. atoms with a specific altloc.

Note for completeness: it is also possible to obtain conformers of an individual residue group. However, conceptually this is best viewed as a shortcut for first obtaining the conformers of a chain, and then finding the residue of interest in each conformer, ignoring all other residues.

• **Residue**: complete residue in a conformer with main conf. atoms and alt. conf. atoms.

Residues are classified as follows:

- *pure main conf.*: all main conf. atoms
- *pure alt. conf.* : all alt. conf. atoms
- *proper alt. conf.*: both main conf. atoms and alt. conf atoms, all alt conf. atoms have a non-blank altloc column
- *improper alt. conf*.: both main conf. atoms and alt. conf atoms, one or more alt conf. atoms have a blank altloc column (as shown in the 1SO7 fragment above).

Errors and warnings

The tools in iotbx.pdb are designed to be as tolerant as reasonably possible when processing input PDB files. E.g., as of Jun 8 2010, all 65802 files in the PDB archive can be processed without generating exceptions. However, to assist users and developers in creating PDB files without ambiguities, the hierarchy object provides methods for flagging likely problems as errors and warnings. It is up to the application how to react to the diagnostics.

The phenix.pdb.hierarchy command can be used to quickly obtain a summary of the hierarchy in a PDB file and some diagnostics. This example command highlights all main features:

```
phenix.pdb.hierarchy pdb1jxw.ent.gz
```

Output:

```
file pdb1jxw.ent.gz
 total number of:
   models:
                 1
   chains:
   alt. conf.: 5
   residues:
   atoms:
                786
   anisou:
                 Ω
 number of atom element+charge types: 5
 histogram of atom element+charge frequency:
   " H " 383
" C " 261
   " O " 77
" N " 59
" S " 6
 residue name classes:
    "common_amino_acid" 48
    "other"
 number of chain ids: 2
 histogram of chain id frequency:
    " " 1
    "A" 1
 number of alt. conf. ids: 3
 histogram of alt. conf. id frequency:
    "A" 2
    "B" 2
    "C" 1
  residue alt. conf. situations:
   pure main conf.:
                      32
   pure alt. conf.:
   proper alt. conf.: 14 improper alt. conf.: 0
 chains with mix of proper and improper alt. conf.: 0
 number of residue names: 16
 histogram of residue name frequency:
   "CYS" 6
   "ALA" 5
    "ILE" 5
    "PRO" 5
    "GLY" 4
    "ASN" 3
    "SER" 3
    "ARG" 2
    "LEU" 2
    "TYR" 2
    "VAL" 2
    "ASP" 1
    "EOH" 1
               other
    "GLU" 1
    "PHE" 1
  ### WARNING: consecutive residue_groups with same resid ###
 number of consecutive residue groups with same resid: 2
   residue group:
            378 N PRO A 22 .*.
      "ATOM
      ... 12 atoms not shown
      "ATOM
              391 HD3APRO A 22 .*.
   next residue group:
      "ATOM 392 CB BSER A 22 .*.
                                          С "
      ... 6 atoms not shown
                                          н "
      "ATOM 399 HB3CSER A 22 .*.
   residue group:
             432 N LEU A 25 .*.
                                          N "
      "ATOM
      ... 18 atoms not shown
      "ATOM 439 CD1CLEU A 25 .*.
                                          C "
   next residue group:
                                          C "
      "ATOM 452 CG1BILE A 25 .*.
      ... 13 atoms not shown
      "ATOM
              466 HD13CILE A 25 .*.
                                          н "
```

The diagnostics are mainly intended for PDB working files, but we have tested the iotbx.pdb module by processing the entire PDB archive. (This took about 3700 CPU seconds, but using 40 CPUs the last job finished after only 277 seconds. Disk and network I/O is rate-limiting in this case, not the performance of the PDB handling library.) The results are summarized in the following table:

```
Total number of .ent files: 65802 (2010 Jun 8)
56873
        WARNING: duplicate chain id
    3
        WARNING: consecutive residue groups with same resid
    2
                 duplicate atom labels
        ERROR:
    1
        ERROR:
                 improper alt. conf.
    0
        ERROR:
                 duplicate model id
        ERROR:
                 residue group with multiple resnames using same altloc
```

About 86% of the PDB entries re-use the same chain id for multiple chains (which are either separated by other chains or TER cards). In many cases, the re-used chain id is the blank character, which is clearly a minor issue. However, we flag this situation to assist people in producing new PDB files with unambiguous chain ids.

The next item in the list, "consecutive residue_groups with same resid" (where resid=resseq+icode), was mentioned before. This situation is best avoided to minimize the chances of mis-interpreting the PDB files.

"duplicate atom labels" pose a serious practical problem (therefore flagged as "ERROR") since it is impossible to uniquely select atoms with duplicate labels, for example via an atom selection syntax as used in many programs (*CNS*, *PyMOL*, *PHENIX*, *VMD*, etc.) or via PDB records in the connectivity annotation section (LINK, SSBOND, CISPEP). Atom serial numbers are not suitable for this purpose since many programs do not preserve them.

Only one PDB entry (1JRT) includes "improper alt. conf." as introduced in the previous section.

There are no PDB entries with two other potential problems diagnosed by the iotbx.pdb module. MODEL ids are unique throughout the PDB archive (as an aside: and all MODEL records have a matching ENDMDL record). Finally, in all 74 files with mixed residue names (i.e. conformers with different sequences), there is exactly one residue name for a given altloc.

Acknowledgments

We gratefully acknowledge the financial support of NIH/NIGMS under grant number P01GM063210. Our work was supported in part by the US Department of Energy under Contract No. DE-AC02-05CH11231.

References

Grosse-Kunstleve, R.W., Zwart, P.H., Afonine, P.V., Ioerger, T.R., Adams, P.D. (2006). Newsletter of the IUCr Commission on Crystallographic Computing, 7, 92-105.